

[KO]MAD2L2 Rabbit mAb


ABZ4631

Product Overview


Name	[KO]MAD2L2 Rabbit mAb
Accession(Primary)	Q9UI95
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (Q9UI95)


Gene Name	MAD2L2
About Gene Symbol	MAD2L2 (Mitotic Arrest Deficient 2 Like 2) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MAD2L2 in HeLa(WT) and MAD2L2-Knockout HeLa(KO) using MAD2L2 antibody.



Western blot analysis of Mad2L2/REV7 in Jurkat lysates using MAD2L2 antibody.



Immunocytochemistry analysis of Mad2L2/REV7 (green) in Hela using Mad2L2/REV7 antibody, and DAPI(blue)