

[KD]NAK/TBK1 Rabbit mAb


ABZ4632

Product Overview


Name	[KD]NAK/TBK1 Rabbit mAb
Accession(Primary)	Q9UHD2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; IP-1:20-1:50

Target information (Q9UHD2)


Gene Name	TBK1
About Gene Symbol	TBK1 (TANK Binding Kinase 1) is a serine/threonine protein kinase involved in NF- κ B signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TBK1 in HeLa(WT) and TBK1-Knockdown HeLa(KD) using TBK1 antibody.



Western blot analysis of TNK1/NAK in Raw264.7, mouse brain lysates using TBK1 antibody.



Western blot analysis of TNK1/NAK in Jurkat, C6, CHO-K1, HeLa lysates using TNK1/NAK antibody.