

[KD]NAK/TBK1 Rabbit mAb


ABZ4634

Product Overview


Name	[KD]NAK/TBK1 Rabbit mAb
Accession(Primary)	Q9UHD2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:200-1:1000; IF: 1:50

Target information (Q9UHD2)


Gene Name	TBK1
About Gene Symbol	TBK1 (TANK Binding Kinase 1) is a serine/threonine protein kinase involved in NF- κ B signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of TBK1 in HeLa(WT) and TBK1-Knockdown HeLa(KD) using TBK1 antibody.



Western blot analysis of TBK1 in MCF-7, SKOV-3 and K562 lysates using TBK1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human testis using NAK/TBK1 (Nterm) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of TBK1 in MCF-7 using NAK/TBK1 (Nterm) antibody.