

[KO]P70 S6 Kinase beta/SRK Rabbit mAb


ABZ4635

Product Overview


Name	[KO]P70 S6 Kinase beta/SRK Rabbit mAb
Accession(Primary)	Q9UBS0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (Q9UBS0)


Gene Name	RPS6KB2
About Gene Symbol	RPS6KB2 (Ribosomal Protein S6 Kinase B2) is a serine/threonine protein kinase involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular m




Western blot analysis of S6K2 in HeLa(WT) and S6K2-Knockout HeLa(KO) using S6K2 antibody.



Immunocytochemistry analysis of S6K2 (green) in 3T3 using S6K2 antibody, and DAPI(blue).



Western blot analysis of P70 S6 Kinase beta in rat Brain, C6, HeLa lysates using P70 S6 Kinase beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using P70 S6 Kinase beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.