

[KO]PAR4 Rabbit mAb


ABZ4638

Product Overview


Name	[KO]PAR4 Rabbit mAb
Accession(Primary)	Q96IZ0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 37 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q96IZ0)


Gene Name	PAWR
About Gene Symbol	PAWR (Pro-Apoptotic WT1 Regulator) is a transcription factor involved in apoptosis. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PAR4 in HeLa(WT) and PAR4-Knockout HeLa(KO) using PAR4 antibody.



Western blot analysis of PAR4 in 3T3 lysates using PAR4 antibody.



Western blot analysis of PAR4 in HeLa, A549, HL-60, U251, U87-MG lysates using PAR4 antibody.