

[KD]DAZAP1 Rabbit mAb


ABZ4639

Product Overview


Name	[KD]DAZAP1 Rabbit mAb
Accession(Primary)	Q96EP5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100; IP-1:50-1:200

Target information (Q96EP5)

Gene Name	DAZAP1
About Gene Symbol	DAZAP1 (DAZ Associated Protein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of DAZAP1 in HeLa(WT) and DAZAP1-Knockdown HeLa(KD) using DAZAP1 antibody.



Western blot analysis of DAZAP1 in K562, C6, 3T3, HeLa lysates using DAZAP1 antibody.