

[KO]Phospho-NDRG1 (Ser330) Rabbit mAb


ABZ4642

Product Overview


Name	[KO]Phospho-NDRG1 (Ser330) Rabbit mAb
Accession(Primary)	Q92597
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:100-1:200; IF-1:50-1:200

Target information (Q92597)


Gene Name	NDRG1
About Gene Symbol	NDRG1 (N-Myc Downstream Regulated 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Phospho-NDRG1 (Ser330) in HeLa(WT) and Phospho-NDRG1 (Ser330)-Knockout HeLa(KO) using Phospho-NDRG1 (Ser330) antibody.




Immunocytochemistry analysis of Phospho-NDRG1 (Ser330) (green) in 293T using Phospho-NDRG1 (Ser330) antibody, and DAPI(blue).



Western blot analysis of Phospho-NDRG1 (Ser330) in mouse kidney lysates using Phospho-NDRG1 (Ser330) antibody.



Western blot analysis of Phospho-NDRG1(Ser330) in Hela, A549, MCF-7, 239T, C2C12 lysates using Phospho-NDRG1(Ser330) antibody



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using NDRG1 (Phospho-Ser330) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.