

[KD]Olig1 Rabbit mAb

ABZ4644


Product Overview

Name	[KD]Olig1 Rabbit mAb
Accession(Primary)	Q8TAK6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100


Target information (Q8TAK6)

Gene Name OLIG1


About Gene Symbol



Western blot analysis of Olig1 in HeLa(WT) and Olig1-Knockdown HeLa(KD) using Olig1 antibody.



Western blot analysis of Olig1 in 3T3 lysates using Olig1 antibody.



Western blot analysis of Olig1 in A549, HL-60, U2OS, C6 lysates using Olig1 antibody.