

[KD]Mitofusin 1 Rabbit mAb

ABZ4645


Product Overview

Name	[KD]Mitofusin 1 Rabbit mAb
Accession(Primary)	Q8R4Z9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50


Target information (Q8R4Z9)

Gene Name Mfn1


About Gene Symbol



Western blot analysis of Mitofusin 1 in HeLa(WT) and Mitofusin 1-Knockdown HeLa(KD) using Mitofusin 1 antibody.



Western blot analysis of Mitofusin 1 in 3T3, mouse brain lysates using Mitofusin 1 antibody.



Western blot analysis of Mitofusin 1 in CHO-K1, C6 lysates using Mitofusin 1 antibody.