

[KD]MLKL Rabbit mAb


ABZ4647

Product Overview


Name	[KD]MLKL Rabbit mAb
Accession(Primary)	Q8NB16
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:50-1:200; ICC/IF-1:100-1:200

Target information (Q8NB16)


Gene Name	MLKL
About Gene Symbol	MLKL (Mixed Lineage Kinase Domain-Like Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of MLKL in HeLa(WT) and MLKL-Knockdown HeLa(KD) using MLKL antibody.




Western blot analysis of MLKL in HUVEC and RAW264.7 lysates using MLKL antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using MLKL antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of MLKL in 3T3, SW480, Ramos and Hela lysates using MLKL antibody.



Western blot analysis of MLKL in HUVEC lysates using MLKL antibody

