

[KD]IGFBP7 Rabbit mAb


ABZ4648

Product Overview


Name	[KD]IGFBP7 Rabbit mAb
Accession(Primary)	Q16270
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:100-1:200; IF-1:50-1:100; FC-1:50-1:100

Target information (Q16270)


Gene Name	IGFBP7
About Gene Symbol	IGFBP7 (Insulin Like Growth Factor Binding Protein 7) is a cytokine and immune signaling protein involved in cell cycle regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of IGFBP7 in HT1080(WT) and IGFBP7-Knockdown HT1080(KD) using IGFBP7 antibody.



Western blot analysis of IGFBP7 in K562, HeLa lysates using IGFBP7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using IGFBP7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.