

[KO]E2F4 Rabbit mAb


ABZ4649

Product Overview


Name	[KO]E2F4 Rabbit mAb
Accession(Primary)	Q16254
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 44 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (Q16254)

Gene Name	E2F4
About Gene Symbol	E2F4 (E2F Transcription Factor 4) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of E2F4 in HeLa(WT) and E2F4-Knockout HeLa(KO) using E2F4 antibody.



Western blot analysis of E2F4 in K562 lysates using E2F4 antibody.