

[KD]G3BP Rabbit mAb


ABZ4655

Product Overview


Name	[KD]G3BP Rabbit mAb
Accession(Primary)	Q13283
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500; FC-1:20-1:50; IP-1:20

Target information (Q13283)


Gene Name	G3BP1
About Gene Symbol	G3BP1 (G3BP Stress Granule Assembly Factor 1) is a GTPase-activating protein involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cel




Western blot analysis of G3BP in HeLa(WT) and G3BP-Knockdown HeLa(KD) using G3BP antibody.



Western blot analysis of G3BP in K562, C6, 3T3, HeLa lysates using G3BP antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using G3BP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of G3BP (green) in Hela using G3BP antibody, and DAPI (blue)