

[KD]G3BP Rabbit mAb


ABZ4656

Product Overview


Name	[KD]G3BP Rabbit mAb
Accession(Primary)	Q13283
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 52 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500; FC-1:20; IP-1:20-1:50

Target information (Q13283)


Gene Name	G3BP1
About Gene Symbol	G3BP1 (G3BP Stress Granule Assembly Factor 1) is a GTPase-activating protein involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cel



Western blot analysis of G3BP in HeLa(WT) and G3BP-Knockdown HeLa(KD) using G3BP antibody.



Western blot analysis of G3BP in K562, C6, HeLa lysates using G3BP antibody.



Immunocytochemistry analysis of G3BP(green) in Hela using G3BP antibody, and DAPI(blue)