

[KO]Phospho-HSF1 (Ser326) Rabbit mAb


ABZ4663

Product Overview


Name	[KO]Phospho-HSF1 (Ser326) Rabbit mAb
Accession(Primary)	Q00613
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 57 kDa; Observed MW: 82 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:50-1:100; IF-1:50-1:200

Target information (Q00613)


Gene Name	HSF1
About Gene Symbol	HSF1 (Heat Shock Transcription Factor 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Phospho-HSF1 (Ser326) in 293T(WT) and Phospho-HSF1 (Ser326)-Knockout 293T(KO) using Phospho-HSF1 (Ser326) antibody.



Western blot analysis of Phospho-HSF1 (Ser326) in K562, HeLa lysates using Phospho-HSF1 (Ser326) antibody.



Immunocytochemistry analysis of HSF1 (Phospho-Ser326)(green) in Hela using HSF1 (Phospho-Ser326) antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using HSF1 (Phospho-Ser326) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.