

[KO]PRAME Rabbit mAb


ABZ4665

Product Overview


Name	[KO]PRAME Rabbit mAb
Accession(Primary)	P78395
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:1000-1:10000; IHC-1:100-1:500; ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:20-1:50

Target information (P78395)


Gene Name	PRAME
About Gene Symbol	



Western blot analysis of PRAME in HCT116(WT) and PRAME-Knockout HCT116(KO) lysates using PRAME antibody.



Immunohistochemistry analysis of paraffin-embedded mouse testis using PRAME antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PRAME in K562 lysates
using PRAME antibody.