

[KD]eIF3B Rabbit mAb


ABZ4667

Product Overview


Name	[KD]eIF3B Rabbit mAb
Accession(Primary)	P55884
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 92 kDa; Observed MW: 117 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:10-1:100

Target information (P55884)


Gene Name	EIF3B
About Gene Symbol	EIF3B (Eukaryotic Translation Initiation Factor 3 Subunit B) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of eIF3B in HeLa(WT) and eIF3B-Knockdown HeLa(KD) using eIF3B antibody.



Western blot analysis of eIF3B in Raji, COS7, 3T3, HeLa lysates using eIF3B antibody.



Immunocytochemistry analysis of eIF3B(green) in HeLa using eIF3B antibody, and DAPI(blue)