

[KD]Presenilin 1 Rabbit mAb



ABZ4674


Product Overview

Name	[KD]Presenilin 1 Rabbit mAb
Accession(Primary)	P49768
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 53 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; IF-1:50-1:200


Target information (P49768)

Gene Name	PSEN1
About Gene Symbol	PSEN1 (Presenilin 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in neurodegenerative disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Western blot analysis of Presenilin 1/PS1 in mouse colon, mouse gastric lysates using Presenilin 1/PS1 antibody.



Western blot analysis of Presenilin 1/PS1 in HeLa, C6, Jurkat lysates using Presenilin 1/PS1 antibody