

[KD]CT A Rabbit mAb


ABZ4676

Product Overview


Name	[KD]CT A Rabbit mAb
Accession(Primary)	P49585
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:20-1:50

Target information (P49585)

Gene Name	PCYT1A
About Gene Symbol	PCYT1A (Phosphate Cytidylyltransferase 1A, Choline) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanism



Western blot analysis of CT A in HeLa(WT) and CT A-Knockdown HeLa(KD) using CT A antibody.



Western blot analysis of CT A in MCF-7, 3T3, HeLa lysates using CT A antibody.