

[KO]PI3 Kinase p110 beta Rabbit mAb


ABZ4680

Product Overview


Name	[KO]PI3 Kinase p110 beta Rabbit mAb
Accession(Primary)	P42338
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 123 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P42338)

Gene Name	PIK3CB
About Gene Symbol	PIK3CB (Phosphatidylinositol-4,5-Bisphosphate 3-Kinase) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely s



Western blot analysis of PI3 Kinase p110 beta in HT1080(WT) and PI3 Kinase p110 beta-Knockout HT1080(KO) using PI3 Kinase p110 beta antibody.



Western blot analysis of PI3 Kinase p110 beta in HCT116, SH-SY5Y lysates using PI3 Kinase p110 beta antibody.