

[KD]FEN1 Rabbit mAb


ABZ4684

Product Overview


Name	[KD]FEN1 Rabbit mAb
Accession(Primary)	P39748
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000

Target information (P39748)

Gene Name	FEN1
About Gene Symbol	FEN1 (Flap Structure-Specific Endonuclease 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of FEN1 in HT1080(WT) and FEN1-Knockdown HT1080(KD) using FEN1 antibody.



Western blot analysis of FEN1 in HeLa, MCF-7, U251, HL-60 lysates using FEN1 antibody.