

[KD]MCM5 Rabbit mAb


ABZ4687

Product Overview


Name	[KD]MCM5 Rabbit mAb
Accession(Primary)	P33992
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 82 kDa; Observed MW: 82 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:2000; IF-1:50-1:200

Target information (P33992)


Gene Name	MCM5
About Gene Symbol	MCM5 (Minichromosome Maintenance Complex Component 5) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechani




Western blot analysis of MCM5 in HeLa(WT) and MCM5-Knockdown HeLa(KD) using MCM5 antibody.



Western blot analysis of MCM5 in Jurkat, C6, HeLa lysates using MCM5 antibody.



Immunocytochemistry analysis of MCM5(green) in Hela using MCM5 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MCM5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.