

[KD]PKC delta Rabbit mAb


ABZ4689

Product Overview


Name	[KD]PKC delta Rabbit mAb
Accession(Primary)	P28867
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 78 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:20-1:500; ICC/IF-1:50-1:100; FC-1:20-1:50

Target information (P28867)


Gene Name	Prkcd
About Gene Symbol	




Western blot analysis of PKC delta in HeLa(WT) and PKC delta-Knockdown HeLa(KD) using PKC delta antibody.



Immunocytochemistry analysis of PKC delta (green) in 3T3 using PKC delta antibody, and DAPI (blue).



Western blot analysis of PKC delta in K562, C6, 3T3 lysates using PKC delta antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using PKC delta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.