

[KD]Phospho-PKC alpha (Ser657) Rabbit mAb


ABZ4691

Product Overview


Name	[KD]Phospho-PKC alpha (Ser657) Rabbit mAb
Accession(Primary)	P17252
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 77 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:1000-1:5000

Target information (P17252)


Gene Name	PRKCA
About Gene Symbol	PRKCA (Protein Kinase C Alpha) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in aging research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Phospho-PKC alpha (Ser657) in HT1080(WT) and Phospho-PKC alpha (Ser657)-Knockdown HT1080(KD) using Phospho-PKC alpha (Ser657) antibody.



Western blot analysis of Phospho-PKC alpha (Ser657) in C6, rat brain lysates using Phospho-PKC alpha (Ser657) antibody.



Western blot analysis of PKC alpha (Phospho-Ser657) in mouse Brain, Ramos lysates using PKC alpha (Phospho-Ser657) antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using Phospho-PKC alpha (Ser657) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.