

[KO]Cyclin E2 Rabbit mAb


ABZ4699

Product Overview


Name	[KO]Cyclin E2 Rabbit mAb
Accession(Primary)	O96020
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:500; ICC/IF-1:50-1:200; FC-1:20; IP-1:20-1:50

Target information (O96020)


Gene Name	CCNE2
About Gene Symbol	CCNE2 (Cyclin E2) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cyclin E2 in HeLa(WT) and Cyclin E2-Knockout HeLa(KO) using Cyclin E2 antibody.



Western blot analysis of Cyclin E2 in K562, HeLa lysates using Cyclin E2 antibody.



Immunocytochemistry analysis of Cyclin
E2(green) in Hela using Cyclin E2 antibody, and
DAPI(blue)