

[KO]FGF2 Rabbit mAb


ABZ4703

Product Overview


Name	[KO]FGF2 Rabbit mAb
Accession(Primary)	P09038
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 31 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:20-1:50

Target information (P09038)

Gene Name	FGF2
About Gene Symbol	FGF2 (Fibroblast Growth Factor 2) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of FGF2 in HeLa(WT) and FGF2-Knockout HeLa(KO) using FGF2 antibody.



Western blot analysis of FGF2 in HeLa lysates using FGF2 antibody.