

[KO]CRABP2 Rabbit mAb


ABZ4704

Product Overview


Name	[KO]CRABP2 Rabbit mAb
Accession(Primary)	P29373
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; FC-1:200-1:500

Target information (P29373)


Gene Name	CRABP2
About Gene Symbol	CRABP2 (Cellular Retinoic Acid Binding Protein 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CRABP2 in HeLa(WT) and CRABP2-Knockout HeLa(KO) using CRABP2 antibody.



Western blot analysis of CRABP2 in HeLa lysates using CRABP2 antibody.



Immunocytochemistry analysis of CRABP2 (green)
in MCF-7 using CRABP2 antibody, and DAPI (blue)