

[KD]Phospho-Presenilin 1 (Ser310) Rabbit mAb


ABZ4707

Product Overview


Name	[KD]Phospho-Presenilin 1 (Ser310) Rabbit mAb
Accession(Primary)	P49768
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 53 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P49768)

Gene Name	PSEN1
About Gene Symbol	PSEN1 (Presenilin 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in neurodegenerative disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-Presenilin 1 (Ser310) in HeLa(WT) and Phospho-Presenilin 1 (Ser310)-Knockdown HeLa(KD) using Phospho-Presenilin 1 (Ser310) antibody.



Western blot analysis of Phospho-Presenilin 1/PS1 (Ser310) in Jurkat, Hela lysates using Phospho-Presenilin 1 (Ser310) antibody.