

[KO]PDGFR beta Rabbit mAb


ABZ4709

Product Overview


Name	[KO]PDGFR beta Rabbit mAb
Accession(Primary)	P09619
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 124 kDa; Observed MW: 190 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:200; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:20

Target information (P09619)


Gene Name	PDGFRB
About Gene Symbol	



Western blot analysis of PDGFR beta in HeLa(WT) and PDGFR beta-Knockout HeLa(KO) using PDGFR beta antibody.



Western blot analysis of PDGFR beta in SH-SY5Y, C6 lysates using PDGFR beta antibody.



Western blot analysis of PDGFR beta in 3T3
lysates using PDGFR beta antibody.