

[KO]Catalase Rabbit mAb


ABZ4711

Product Overview


Name	[KO]Catalase Rabbit mAb
Accession(Primary)	P24270
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50-1:100; FC-1:50-1:100

Target information (P24270)


Gene Name	Cat
About Gene Symbol	



Western blot analysis of Catalase in HeLa(WT) and Catalase-Knockout HeLa(KO) using Catalase antibody.



Western blot analysis of Catalase in K562, rat Brain, C6, 3T3 lysates using Catalase antibody.



Immunocytochemistry analysis of Catalase
(green) in Raw264.7 using Catalase antibody, and
DAPI(blue)