

[KD]Phospho-CDC6 (Ser106) Rabbit mAb


ABZ4712

Product Overview


Name	[KD]Phospho-CDC6 (Ser106) Rabbit mAb
Accession(Primary)	Q99741
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 63 kDa; Observed MW: 63 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q99741)

Gene Name	CDC6
About Gene Symbol	CDC6 (Cell Division Cycle 6) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-CDC6 (Ser106) in HeLa(WT) and Phospho-CDC6 (Ser106)-Knockdown HeLa(KD) using Phospho-CDC6 (Ser106) antibody.



Western blot analysis of Phospho-Cdc6 (Ser106) in Jurkat lysates using Phospho-CDC6 (Ser106) antibody.