

[KD]Creatine Kinase B type Rabbit mAb


ABZ4714

Product Overview


Name	[KD]Creatine Kinase B type Rabbit mAb
Accession(Primary)	P12277
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:1000; ICC/IF-1:50-1:500; FC-1:10-1:100

Target information (P12277)


Gene Name	CKB
About Gene Symbol	CKB (Creatine Kinase B) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Creatine Kinase B type in HeLa(WT) and Creatine Kinase B type-Knockdown HeLa(KD) using Creatine Kinase B type antibody.



Western blot analysis of Creatine Kinase B type in mouse brain lysates using Creatine Kinase B type antibody.



Western blot analysis of Creatine kinase B type in K562, rat Brain, C6, HeLa lysates using Creatine kinase B type antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Creatine kinase B type antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.