

[KD]Cyclin E1 Rabbit mAb


ABZ4715

Product Overview


Name	[KD]Cyclin E1 Rabbit mAb
Accession(Primary)	P24864
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:500-1:1000; FC-1:50-1:100

Target information (P24864)

Gene Name	CCNE1
About Gene Symbol	CCNE1 (Cyclin E1) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cyclin E1 in HeLa(WT) and Cyclin E1-Knockdown HeLa(KD) using Cyclin E1 antibody.



Western blot analysis of Cyclin E1 in K562, HeLa lysates using Cyclin E1 antibody.