

[KD]DARPP32 Rabbit mAb


ABZ4716

Product Overview


Name	[KD]DARPP32 Rabbit mAb
Accession(Primary)	Q9UD71
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:50-1:100; IP-1:20

Target information (Q9UD71)


Gene Name	PPP1R1B
About Gene Symbol	




Western blot analysis of DARPP32 in HeLa(WT) and DARPP32-Knockdown HeLa(KD) using DARPP32 antibody.



Western blot analysis of DARPP32 in mouse brain lysates using DARPP32 antibody.



Western blot analysis of DARPP32 in rat Brain lysates using DARPP32 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using DARPP32 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.