

[KD]DKC1 Rabbit mAb


ABZ4718

Product Overview


Name	[KD]DKC1 Rabbit mAb
Accession(Primary)	O60832
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200

Target information (O60832)


Gene Name	DKC1
About Gene Symbol	




Western blot analysis of DKC1 in HeLa(WT) and DKC1-Knockdown HeLa(KD) using DKC1 antibody.




Immunocytochemistry analysis of DKC1 (green) in K562 using DKC1 antibody, and DAPI(blue).





Western blot analysis of DKC1 in K562, C6, 3T3, HeLa lysates using DKC1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using DKC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.