

[KO]LAT Rabbit mAb


ABZ4721

Product Overview


Name	[KO]LAT Rabbit mAb
Accession(Primary)	O43561
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 28 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; FC-1:100

Target information (O43561)


Gene Name	LAT
About Gene Symbol	LAT (Linker For Activation Of T Cells) is a signaling adaptor protein involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of LAT in HT1080(WT) and LAT-Knockout HT1080(KO) using LAT antibody.



Western blot analysis of LAT in U2OS lysates using LAT antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using LAT antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.