

[KO]PDPK1 Rabbit mAb


ABZ4722

Product Overview


Name	[KO]PDPK1 Rabbit mAb
Accession(Primary)	O15530
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 63 kDa; Observed MW: 63 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (O15530)

Gene Name	PDPK1
About Gene Symbol	PDPK1 (3-Phosphoinositide Dependent Protein Kinase 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely stu



Western blot analysis of PDPK1 in HeLa(WT) and PDPK1-Knockout HeLa(KO) using PDPK1 antibody.



Western blot analysis of PDPK1 in HeLa lysates using PDPK1 antibody.