

[KD]Phospho-PKC gamma (Thr514) Rabbit mAb


ABZ4723

Product Overview


Name	[KD]Phospho-PKC gamma (Thr514) Rabbit mAb
Accession(Primary)	P05129
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 78 kDa; Observed MW: 78-85 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:20

Target information (P05129)


Gene Name	PRKCG
About Gene Symbol	PRKCG (Protein Kinase C Gamma) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Phospho-PKC gamma (Thr514) in HT1080(WT) and Phospho-PKC gamma (Thr514)-Knockdown HT1080(KD) using Phospho-PKC gamma (Thr514) antibody.



Western blot analysis of Phospho-PKC gamma (Thr514) in mouse brain lysates using Phospho-PKC gamma (Thr514) antibody.



Western blot analysis of Phospho-PKC gamma (Thr514) in C6, HeLa lysates using Phospho-PKC gamma (Thr514) antibody.



Immunohistochemistry analysis of paraffin-embedded Human Brain using PKC (Phospho-Thr514) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.