

[KO]Phospho-PKR (Thr451) Rabbit mAb


ABZ4724

Product Overview


Name	[KO]Phospho-PKR (Thr451) Rabbit mAb
Accession(Primary)	P19525
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 62 kDa; Observed MW: 74 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P19525)

Gene Name	EIF2AK2
About Gene Symbol	EIF2AK2 (eIF-2A protein kinase 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-PKR (Thr451) in HeLa(WT) and Phospho-PKR (Thr451)-Knockout HeLa(KO) using Phospho-PKR (Thr451) antibody.



Western blot analysis of PKR (Phospho-Thr451) in HeLa, HeLa+CalyculinA lysates using Phospho-PKR (Thr451) antibody.