

[KO]Placental Alkaline Phosphatase Rabbit mAb


ABZ4725

Product Overview


Name	[KO]Placental Alkaline Phosphatase Rabbit mAb
Accession(Primary)	P05187
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 58 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:200-1:500

Target information (P05187)


Gene Name	ALPP
About Gene Symbol	ALPP (Alkaline Phosphatase, Placental) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Placental Alkaline Phosphatase in HeLa(WT) and Placental Alkaline Phosphatase-Knockout HeLa(KO) using Placental Alkaline Phosphatase antibody.



Western blot analysis of Placental alkaline phosphatase (Phospho- LAP) in HeLa lysates using Placental Alkaline Phosphatase antibody.



Western blot analysis of Placental alkaline phosphatase (PLAP) in HeLa lysates using Placental alkaline phosphatase (PLAP) antibody.