

[KD]Rab5 Rabbit mAb


ABZ4726

Product Overview


Name	[KD]Rab5 Rabbit mAb
Accession(Primary)	P20339
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:500-1:1000; FC-1:20-1:50; IP-1:20-1:50

Target information (P20339)


Gene Name	RAB5A
About Gene Symbol	RAB5A (RAB5A, Member RAS Oncogene Family) is a small GTPase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Rab5 in HeLa(WT) and Rab5-Knockdown HeLa(KD) using Rab5 antibody.




Immunocytochemistry analysis of Rab5 (green) in HT-1080 using Rab5 antibody, and DAPI(blue).





Western blot analysis of Rab5 in K562, C6, 3T3, HeLa lysates using Rab5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using Rab5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.