

## [KD]Src Rabbit mAb


ABZ4728

### Product Overview


|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | [KD]Src Rabbit mAb                         |
| <b>Accession(Primary)</b>   | P12931                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 60 kDa; Observed MW: 60 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:10000; IHC-1:50-1:2000         |

### Target information (P12931)


|                          |                                                                                                                                                                                                                                                                                     |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | SRC                                                                                                                                                                                                                                                                                 |
| <b>About Gene Symbol</b> | <p>SRC (SRC Proto-Oncogene, Non-Receptor Tyrosine Kinase) is a receptor tyrosine kinase involved in receptor-mediated signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.</p> |




Western blot analysis of Src in HT1080(WT) and Src-Knockdown HT1080(KD) using Src antibody.



Western blot analysis of Src in SH-SY5Y, U87-MG lysates using Src antibody.



Western blot analysis of Src in C6, 3T3 lysates using Src antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Src antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.