

[KO]Annexin 4 Rabbit mAb


ABZ4732

Product Overview


Name	[KO]Annexin 4 Rabbit mAb
Accession(Primary)	P09525
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:200-1:500

Target information (P09525)


Gene Name	ANXA4
About Gene Symbol	ANXA4 (Annexin A4) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Annexin 4 in HeLa(WT) and Annexin 4-Knockout HeLa(KO) using Annexin 4 antibody.




Immunocytochemistry analysis of Annexin 4 (green) in HeLa using Annexin 4 antibody, and DAPI(blue).



Western blot analysis of Annexin4/ANXA4 in mouse kidney, mouse lung lysates using Annexin4/ANXA4 antibody.



Western blot analysis of Annexin A4 in HeLa, CHO-K1, C6 lysates using Annexin A4 antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Annexin4/ANXA4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.