

[KD]Bcl 10 Rabbit mAb


ABZ4734

Product Overview


Name	[KD]Bcl 10 Rabbit mAb
Accession(Primary)	O95999
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200

Target information (O95999)

Gene Name	BCL10
About Gene Symbol	BCL10 (BCL10 Immune Signaling Adaptor) is a cysteine protease (caspase) involved in NF- κ B signaling and apoptosis. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Bcl10 in HeLa(WT) and Bcl10-Knockdown HeLa(KD) using Bcl10 antibody.



Western blot analysis of Bcl10 in A549, MCF-7, 293T lysates using Bcl10 antibody.