

[KO]CDC123 Rabbit mAb


ABZ4735

Product Overview


Name	[KO]CDC123 Rabbit mAb
Accession(Primary)	O75794
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100

Target information (O75794)

Gene Name	CDC123
About Gene Symbol	CDC123 (Cell Division Cycle 123) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CDC123 in HeLa(WT) and CDC123-Knockout HeLa(KO) using CDC123 antibody.



Western blot analysis of CDC123 in HeLa, A549, HL-60 lysates using CDC123 antibody.