

[KD]Cyclin B1 Rabbit mAb


ABZ4741

Product Overview


Name	[KD]Cyclin B1 Rabbit mAb
Accession(Primary)	P24860
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 48 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; FC-1:100-1:200

Target information (P24860)


Gene Name	Ccnb1
About Gene Symbol	



Western blot analysis of Cyclin B1 in HT1080(WT) and Cyclin B1-Knockdown HT1080(KD) using Cyclin B1 antibody.



Western blot analysis of Cyclin B1 in Jurkat lysates using Cyclin B1 antibody.



Immunocytochemistry analysis of Cyclin
B1(green) in Hela using Cyclin B1 antibody, and
DAPI(blue)