

PI3 Kinase p85 beta (8D9) Mouse mAb

ABZ4747


Product Overview

Name	PI3 Kinase p85 beta (8D9) Mouse mAb
Accession(Primary)	O00459
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 82 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O00459)

Gene Name PIK3R2

About Gene Symbol



Western blot analysis of PI3 Kinase p85 beta in C6, K562, MCF-7, U87-MG, HCT116, 3T3 and HeLa lysates using PI3 Kinase p85 beta antibody.

