

PRMT1 (4E9) Mouse mAb


ABZ4748

Product Overview

Name	PRMT1 (4E9) Mouse mAb
Accession(Primary)	Q99873
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q99873)

Gene Name	PRMT1
About Gene Symbol	PRMT1 (Protein Arginine Methyltransferase 1) is a methyltransferase involved in JAK-STAT/cytokine signaling and chromatin remodeling and epigenetics. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PRMT1 in HeLa, A549 and SW480 lysates using PRMT1 antibody.