

IGF2BP1 Rabbit mAb

ABZ4749


Product Overview

Name	IGF2BP1 Rabbit mAb
Accession(Primary)	Q9NZI8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 63 kDa; Observed MW: 63 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:50-1:1000; IP-1:20-1:50


Target information (Q9NZI8)

Gene Name IGF2BP1

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded mouse colon using IGF2BP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of IGF2BP1 in K562, 3T3
lysates using IGF2BP1 antibody.