

Cleaved-Notch1 Rabbit pAb


ABZ4751

Product Overview

Name	Cleaved-Notch1 Rabbit pAb
Accession(Primary)	P46531
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 273 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300; IF-1:50-1:200

Target information (P46531)

Gene Name	NOTCH1
About Gene Symbol	NOTCH1 (Notch Receptor 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cleaved-Notch1 in various lysates using Cleaved-Notch1 antibody.